



polA Antibody, Biotin conjugated

Product Code	CSB-PA16667D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19821
Immunogen	Recombinant Thermus aquaticus DNA polymerase I, thermostable protein (310-832AA)
Raised In	Rabbit
Species Reactivity	Thermus aquaticus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA polymerase I, thermostable (EC 2.7.7.7) (Taq polymerase 1), polA, pol1
Immunogen Species	Thermus aquaticus
Research Area	Others
Target Names	polA