



Ngf Antibody, HRP conjugated

Product Code	CSB-PA16859B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25427
Immunogen	Recombinant Rat Beta-nerve growth factor protein (124-241AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-nerve growth factor (Beta-NGF), Ngf, Ngfb
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Neuroscience
Target Names	Ngf