



IFNG Antibody, HRP conjugated

Product Code	CSB-PA15594B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49708
Immunogen	Recombinant Chicken Interferon gamma protein (21-164AA)
Raised In	Rabbit
Species Reactivity	Chicken
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon gamma (IFN-gamma), IFNG
Immunogen Species	Chicken
Research Area	Immunology
Target Names	IFNG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.