



# cagA Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA14367C0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P55980
<b>Immunogen</b>	Recombinant Helicobacter pylori Cytotoxicity-associated immunodominant antigen protein (1-294AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Helicobacter pylori
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Cytotoxicity-associated immunodominant antigen (120 kDa protein) (CAG pathogenicity island protein 26), cagA, cag26 cai
<b>Immunogen Species</b>	Helicobacter pylori (strain ATCC 700392 / 26695) (Campylobacter pylori)
<b>Research Area</b>	Others
<b>Target Names</b>	cagA
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.