




Tp53 Antibody

Product Code	CSB-PA15509A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10361
Immunogen	Recombinant Rat Cellular tumor antigen p53 protein (1-391AA)
Raised In	Rabbit
Species Reactivity	Rat, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cellular tumor antigen p53 (Tumor suppressor p53), Tp53, P53
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cell Biology
Target Names	Tp53

Image



Western Blot

Positive WB detected in: Mouse kidney tissue

All lanes: Tp53 antibody at 2.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 44 kDa

Observed band size: 44 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.