



# Recombinant Coprinopsis cinerea Flap endonuclease 1 (FEN1)

<b>Product Code</b>	CSB-BP008585DWQ
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A8NQC2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Coprinopsis cinerea (strain Okayama-7 / 130 / ATCC MYA-4618 / FGSC 9003) (Inky cap fungus) (Hormographiella aspergillata)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGIKGLTGLL NEHAPNSIKE HDIKTLFGRK VAIDASMSIY QFLIAVRQRD GEMLTNDAGE TTSHLMGFFY RTIRIVENGI KPAYVFDGKP PELKKGVLK RFEKREEAKE EGEEAKEIGT AEDVDRFSRR TVKVTKQHNE ECQKLLRLMG VPCVIAPSEA EAQCAELARG GKVYAAGSED MDTLTFNAPI LFRHLTFSEA KKQPISEINL EAALKGLDMD MSQFVDLCIL LGCDYLEPIK GVGPKSALKL IREFGGLKEV VEHLREKAAA RKAEAEDEE EEAEEPAPTS DVEMPDDDEDG EKDSDEEEA ERRKKAEEAAK KKAQEKAKS AKKKGKKGKGIQIPDEWPW EEAKQIFLKP DVIPADQVEL EWKNPDVEGL VQFLVTEKGF SEERVRKGAE KLTFLNAKQ QGRLDGFFTV KPKQAPPAA KGKGAASKGT KRKGEDKAEG SGKKAKKK
<b>Source</b>	Baculovirus
<b>Target Names</b>	FEN1
<b>Protein Names</b>	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
<b>Expression Region</b>	1-458
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.