



# Recombinant Coprinopsis cinerea Flap endonuclease 1 (FEN1)

<b>Product Code</b>	CSB-EP008585DWQ-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A8NQC2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Coprinopsis cinerea (strain Okayama-7 / 130 / ATCC MYA-4618 / FGSC 9003) (Inky cap fungus) (Hormographiella aspergillata)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>           MGIKGLTGLL NEHAPNSIKE HDIKTLFGRK VAIDASMSIY QFLIAVRQRD            GEMLTNDAGE TTSHLMGFFY RTIRIVENGI KPAYVFDGKP PELKKGVLK            RFEKREEAKE EGEEAKEIGT AEDVDRFSRR TVKVTKQHNE ECQKLLRLMG            VPCVIAPSEA EAQCAELARG GKVYAAGSED MDLTLFNAPI LFRHLTFSEA            KKQPISEINL EAALKGLDMD MSQFVDLCIL LGCDYLEPIK GVGPKSALKL            IREFGGLKEV VEHLREKAAA RKAEAEDEE EEAEEPAPTS DVEMPDDDEDG            EKSDSDEEEA ERRKKAEEAAK KKAQEKAKS AKKGKGGKGGK GIQIPDEWPW            EEAKQIFLKP DVIPADQVEL EWKNPDVEGL VQFLVTEKGF SEERVKGAEG            KLTKFLNAKQ QGRLDGFFTV KPKQAPPPAA KGKGAASKGT KRKGEDKAEG            SGKKAKKK         </p>
<b>Source</b>	E.coli
<b>Target Names</b>	FEN1
<b>Protein Names</b>	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
<b>Expression Region</b>	1-458
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.