



Recombinant *Coprinopsis cinerea* Flap endonuclease 1 (FEN1)

Product Code	CSB-MP008585DWQ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A8NQC2
Product Type	Recombinant Protein
Immunogen Species	<i>Coprinopsis cinerea</i> (strain Okayama-7 / 130 / ATCC MYA-4618 / FGSC 9003) (Inky cap fungus) (<i>Hormographiella aspergillata</i>)
Purity	>85% (SDS-PAGE)
Sequence	<p> MGIKGLTGLL NEHAPNSIKE HDIKTLFGRK VAIDASMSIY QFLIAVRQRD GEMLTNDAGE TTSHLMGFFY RTIRIVENGI KPAYVFDGKP PELKKGVLSK RFEKREEAKE EGEEAKEIGT AEDVDRFSRR TVKVTKQHNE ECQKLLRLMG VPCVIAPSEA EAQCAELARG GKVYAAGSED MDLTLFNAPI LFRHLTFSEA KKQPISEINL EAALKGLDMD MSQFVDLCIL LGCDYLEPIK GVGPKSALKL IREFGGLKEV VEHLREKAAA RKAEAEDEE EEAEEPAPTS DVEMPDDDEDG EKSDSDEEEA ERRKKAEEAAK KKAQEKAKS AKKGKGGKGG GIQIPDEWPW EEAKQIFLKP DVIPADQVEL EWKNPDVEGL VQFLVTEKGF SEERVKGAEG KLTKFLNAKQ QGRLDGFFTV KPKQAPPAA KGKGAASKGT KRKGEDKAEG SGKKAKKK </p>
Source	Mammalian cell
Target Names	FEN1
Protein Names	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
Expression Region	1-458
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.