





BMP6 Antibody

Product Code	CSB-PA10709A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22004
Immunogen	Recombinant Human Bone morphogenetic protein 6 protein (382-513AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bone morphogenetic protein 6 (BMP-6) (VG-1-related protein) (VG-1-R) (VGR-1), BMP6, VGR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BMP6

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA10709A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA10709A0Rb at dilution of 1:100

