



# Recombinant Danio rerio Cyclin-dependent kinase 20 (cdk20)

<b>Product Code</b>	CSB-MP004849DIL
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A8WIP6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDQYSILGRI GEGAHGIVFK AKHIETGETV ALKKVALRRL EDGIPNQALR EIKALQEIED NQYVVKLKDV FPHGTGFVLV FEYMLSDLSE VIRNSQRPLT ASQVKSMMMM LLKGVAFCHE NSIMHRDLKP ANLLISSTGH LKIADFGLAR LFSNEGDRLY SHQVATRWYR APELLYGARK YDEGVDLWAV GCIFGELLNN SPLFPGENDI EQLCCVLRVL GTPNQKVWPE ITELPDYNKI TFKENPPIPL EEIVPDTSPQ AVDLLKKFLV YPSKQRISAR QALLHPYFFT DPLPAHHSEL PIPQRGGKHS RQRMQPPHEF TVDRPLHESV VDPSLIQKHA MSCS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	cdk20
<b>Protein Names</b>	Recommended name: Cyclin-dependent kinase 20 EC= 2.7.11.22 Alternative name(s): Cell cycle-related kinase Cell division protein kinase 20
<b>Expression Region</b>	1-344
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.