



mpt63 Antibody, Biotin conjugated

| | |
|---------------------|--|
| Product Code | CSB-PA14907D0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P9WIP1 |
| Immunogen | Recombinant Mycobacterium tuberculosis Immunogenic protein MPT63 protein (30-159AA) |
| Raised In | Rabbit |
| Species Reactivity | Mycobacterium tuberculosis |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Immunogenic protein MPT63 (16 kDa immunoprotective extracellular protein) (Antigen MPT63), mpt63 |
| Immunogen Species | Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv) |
| Research Area | Others |
| Target Names | mpt63 |