



TNNT2 Antibody, Biotin conjugated

Product Code	CSB-PA10047D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45379
Immunogen	Recombinant Human Troponin T, cardiac muscle protein (138-298AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Troponin T, cardiac muscle (TnTc) (Cardiac muscle troponin T) (cTnT), TNNT2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TNNT2