



EGF Antibody, Biotin conjugated

Product Code	CSB-PA06244D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01133
Immunogen	Recombinant Human Pro-epidermal growth factor protein (977-1023AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pro-epidermal growth factor (EGF) [Cleaved into: Epidermal growth factor (Urogastrone)], EGF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EGF