



poxB Antibody

Product Code	CSB-PA17967A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07003
Immunogen	Recombinant Escherichia coli Pyruvate dehydrogenase [ubiquinone] protein (1-572AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyruvate dehydrogenase [ubiquinone] (EC 1.2.5.1) (Pyruvate oxidase) (POX) [Cleaved into: Alpha-peptide], poxB
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	poxB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.