




IL18BP Antibody

Product Code	CSB-PA16739A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95998
Immunogen	Recombinant Human Interleukin-18-binding protein (31-183AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-18-binding protein (IL-18BP) (Tadekinig-alfa), IL18BP
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL18BP

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA16739A0Rb at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.