




Hp Antibody

Product Code	CSB-PA17069A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06866
Immunogen	Recombinant Rat Haptoglobin protein (19-346AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Haptoglobin (Liver regeneration-related protein LRRG173) (Zonulin) [Cleaved into: Haptoglobin alpha chain; Haptoglobin beta chain], Hp
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Hp

Image



Western blot


All lanes: Hp antibody at 2µg/ml + Rat live tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 39, 43 kDa

Observed band size: 39 kDa



Immunohistochemistry of paraffin-embedded mouse liver tissue using CSB-PA17069A0Rb at dilution of 1:50

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.