



# L1 Antibody

Product Code	CSB-PA08359A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06794
Immunogen	Recombinant HPV18 Major capsid protein L1 protein (1-568AA)
Raised In	Rabbit
Species Reactivity	Human papillomavirus type 18
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Human papillomavirus type 18
Research Area	Microbiology
Target Names	L1