



ENPEP Antibody, Biotin conjugated

Product Code	CSB-PA14759D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07075
Immunogen	Recombinant Human Glutamyl aminopeptidase protein (719-949AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamyl aminopeptidase (EAP) (EC 3.4.11.7) (Aminopeptidase A) (AP-A) (Differentiation antigen gp160) (CD antigen CD249), ENPEP
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ENPEP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.