



AGER Antibody, HRP conjugated

Product Code	CSB-PA14209B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15109
Immunogen	Recombinant Human Advanced glycosylation end product-specific receptor protein (23-342AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Advanced glycosylation end product-specific receptor (Receptor for advanced glycosylation end products), AGER, RAGE
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	AGER
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.