



ALPI Antibody, HRP conjugated

Product Code	CSB-PA11989B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09923
Immunogen	Recombinant Human Intestinal-type alkaline phosphatase protein (26-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intestinal-type alkaline phosphatase (IAP) (Intestinal alkaline phosphatase) (EC 3.1.3.1), ALPI
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ALPI
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.