



RELT Antibody, Biotin conjugated

Product Code	CSB-PA08179D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969Z4
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 19L protein (26-153AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 19L (Receptor expressed in lymphoid tissues), RELT, TNFRSF19L
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RELT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.