



Fgf2 Antibody, HRP conjugated

Product Code	CSB-PA06259B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15655
Immunogen	Recombinant Mouse Fibroblast growth factor 2 protein (10-154AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor 2 (FGF-2) (Basic fibroblast growth factor) (bFGF) (Heparin-binding growth factor 2) (HBGF-2), Fgf2, Fgf-2
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Fgf2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.