



# speA Antibody, HRP conjugated

Product Code	CSB-PA14507B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21170
Immunogen	Recombinant Escherichia coli Biosynthetic arginine decarboxylase protein (3-656AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Biosynthetic arginine decarboxylase (ADC) (EC 4.1.1.19), speA
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	speA