



ORM2 Antibody, HRP conjugated

Product Code	CSB-PA14549B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19652
Immunogen	Recombinant Human Alpha-1-acid glycoprotein 2 protein (19-201AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-1-acid glycoprotein 2 (AGP 2) (Orosomuroid-2) (OMD 2), ORM2, AGP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ORM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.