



# IL1RN Antibody, Biotin conjugated

|                            |                                                                                                                              |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA06677D0Rb                                                                                                              |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                |
| <b>Uniprot No.</b>         | P18510                                                                                                                       |
| <b>Immunogen</b>           | Recombinant Human Interleukin-1 receptor antagonist protein (26-177AA)                                                       |
| <b>Raised In</b>           | Rabbit                                                                                                                       |
| <b>Species Reactivity</b>  | Human                                                                                                                        |
| <b>Tested Applications</b> | ELISA                                                                                                                        |
| <b>Form</b>                | Liquid                                                                                                                       |
| <b>Conjugate</b>           | Biotin                                                                                                                       |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                                             |
| <b>Purification Method</b> | >95%, Protein G purified                                                                                                     |
| <b>Isotype</b>             | IgG                                                                                                                          |
| <b>Clonality</b>           | Polyclonal                                                                                                                   |
| <b>Alias</b>               | Interleukin-1 receptor antagonist protein (IL-1RN) (IL-1ra) (IRAP) (ICIL-1RA) (IL1 inhibitor) (Anakinra), IL1RN, IL1F3 IL1RA |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                         |
| <b>Research Area</b>       | Immunology                                                                                                                   |
| <b>Target Names</b>        | IL1RN                                                                                                                        |