



IAPP Antibody, Biotin conjugated

Product Code	CSB-PA12224D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10997
Immunogen	Recombinant Human Islet amyloid polypeptide protein (34-70AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Islet amyloid polypeptide (Amylin) (Diabetes-associated peptide) (DAP) (Insulinoma amyloid peptide), IAPP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IAPP