



IAPP Antibody, Biotin conjugated

Product Code	CSB-PA12224D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10997
Immunogen	Recombinant Human Islet amyloid polypeptide protein (34-70AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Islet amyloid polypeptide (Amylin) (Diabetes-associated peptide) (DAP) (Insulinoma amyloid peptide), IAPP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IAPP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.