



CALCA Antibody, Biotin conjugated

Product Code	CSB-PA12134D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06881
Immunogen	Recombinant Human Calcitonin gene-related peptide 1 protein (83-118AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcitonin gene-related peptide 1 (Alpha-type CGRP) (Calcitonin gene-related peptide I) (CGRP-I), CALCA, CALC1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CALCA