



Recombinant Yersinia pestis bv. Antiqua Phosphatidylserine decarboxylase proenzyme (psd)

Product Code	CSB-MP018033YAM
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A9QYP0
Product Type	Recombinant Protein
Immunogen Species	Yersinia pestis bv. Antiqua (strain Angola)
Purity	>85% (SDS-PAGE)
Sequence	MLDSIKIKLQ YLLPKQGLTQ LAGWGANKQG GWLTQLVIKA FARYYKVDMK EAQDPEFSAY RTFNEFFVRP LRAQVVRPVA EENLLAQPAD GAISQLGAIR EGQILQAKGH NYSLEALLAG NYLLAAEFQN GQFVTTYLAP RDYHRVHMPC DGVLREMIYV PGDLFSVNPL TAANVPNLFA RNERVICIFD TAFGPMAQIL VGATIVGSIE TVWAGTITPP REGVIRRWTY PQAGCEGAIT LEKGQEMGRF KLG
Source	Mammalian cell
Target Names	psd
Protein Names	Recommended name: Phosphatidylserine decarboxylase proenzyme EC= 4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.