



# KRAS Antibody, FITC conjugated

|                            |                                                                                                                            |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA13969C0Rb                                                                                                            |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                              |
| <b>Uniprot No.</b>         | P01116                                                                                                                     |
| <b>Immunogen</b>           | Recombinant Human GTPase KRas protein (2-168AA)                                                                            |
| <b>Raised In</b>           | Rabbit                                                                                                                     |
| <b>Species Reactivity</b>  | Human                                                                                                                      |
| <b>Form</b>                | Liquid                                                                                                                     |
| <b>Conjugate</b>           | FITC                                                                                                                       |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                                           |
| <b>Purification Method</b> | >95%, Protein G purified                                                                                                   |
| <b>Isotype</b>             | IgG                                                                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                                                                 |
| <b>Alias</b>               | GTPase KRas (K-Ras 2) (Ki-Ras) (c-K-ras) (c-Ki-ras) [Cleaved into: GTPase KRas, N-terminally processed], KRAS, KRAS2 RASK2 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                       |
| <b>Research Area</b>       | Signal Transduction                                                                                                        |
| <b>Target Names</b>        | KRAS                                                                                                                       |