



CXCL10 Antibody

Product Code	CSB-PA06107A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02778
Immunogen	Recombinant Human C-X-C motif chemokine 10 protein (22-98AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-X-C motif chemokine 10 (10 kDa interferon gamma-induced protein) (Gamma-IP10) (IP-10) (Small-inducible cytokine B10) [Cleaved into: CXCL10(1-73)], CXCL10, INP10 SCYB10
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CXCL10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.