



LOX Antibody, HRP conjugated

Product Code	CSB-PA16909B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28300
Immunogen	Recombinant Human Protein-lysine 6-oxidase protein (174-417AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein-lysine 6-oxidase (EC 1.4.3.13) (Lysyl oxidase), LOX
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LOX