



# Recombinant Human 3-oxo-5-beta-steroid 4-dehydrogenase (AKR1D1)

<b>Product Code</b>	CSB-BP001548HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P51857
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDLSAASHRI PLSDGNSIPI IGLGTYSEPK STPKGACATS VKVAIDTGYR HIDGAYIQN EHEVGEAIRE KIAEGKVRRE DIFYCGKLWA TNHVP MVRP TLERTLRVLQ LDYVDLYIIE VPMFAKPGDE IYPRDENGKW LYHKS NL CAT WEAMEACKDA GLVKSLGVS N FNRRQLELIL NKPGLKHKPV SNQVECHPYF TQP KLLKFCQ QHDIVITAYS PLGTSRNP IW VNVSSPPLLK DALLNSLGKR YNK TAAQIVL RFNIQRGVVV IPKSFNLERI KENFQIFDFS LTEEEMKDIE ALNKNVRFVE LLMWRDHPEY PFHDEY
<b>Source</b>	Baculovirus
<b>Target Names</b>	AKR1D1
<b>Protein Names</b>	Recommended name: 3-oxo-5-beta-steroid 4-dehydrogenase EC= 1.3.1.3 Alternative name(s): Aldo-keto reductase family 1 member D1 Delta(4)-3-ketosteroid 5-beta-reductase Delta(4)-3-oxosteroid 5-beta-reductase
<b>Expression Region</b>	1-326
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	The enzyme encoded by this gene is responsible for the catalysis of the 5-beta-reduction of bile acid intermediates and steroid hormones carrying a delta(4)-3-one structure. Deficiency of this enzyme may contribute to hepatic dysfunction.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.