



ANXA3 Antibody, HRP conjugated

Product Code	CSB-PA17149B0Rb
Abbreviation	Annexin A3 protein
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12429
Immunogen	Recombinant Human Annexin A3 protein (5-322AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Inhibitor of phospholipase A2, also possesses anti-coagulant properties. Also cleaves the cyclic bond of inositol 1,2-cyclic phosphate to form inositol 1-phosphate.
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Annexin A3 (35-alpha calcimedlin) (Annexin III) (Annexin-3) (Inositol 1,2-cyclic phosphate 2-phosphohydrolase) (Lipocortin III) (Placental anticoagulant protein III) (PAP-III), ANXA3, ANX3
Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ANXA3