



ALDOA Antibody, Biotin conjugated

Product Code	CSB-PA00015D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04075
Immunogen	Recombinant Human Fructose-bisphosphate aldolase A protein (2-364AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fructose-bisphosphate aldolase A (EC 4.1.2.13) (Lung cancer antigen NY-LU-1) (Muscle-type aldolase), ALDOA, ALDA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ALDOA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.