



# Recombinant Human Heterogeneous nuclear ribonucleoprotein M

<b>Product Code</b>	CSB-EP010613HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P52272
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>AAGVEAAAE VAATEIKMEE ESGAPGVPSG NGAPGPKGEG ERPAQNEKRR  EKNIKRGGNR FEPYANPTKR YRAFITNIPF DVKWQSLKDL VKEKVGEVTY  VELLMDAEGK SRGCAVVEFK MEESMKKAAE VLNKHSLSGR PLKVKEDPDG  EHARRAMQKV MATTGGMGMG PGGPGMITIP PSILNNPNIP NEIIHALQAG  RLGSTVFVAN LDYKVGWKKL KEVFSMAGVV VRADILEDKD GKSRGIGTVT  FEQSIEAVQA ISMFNGQLLF DRPMHVKMDE RALPKGDFFP PERPQQLPHG  LGGIGMGLGP GGQPIDANHL NKGIGMGNIG PAGMGMEGIG FGINKMGGME  GPFGGGMENM GRFGSGMNMG RINEILSNAL KRGEIIAKQG  GGGGGGSVPG IERMGPIDR LGGAGMERMG AGLGHGMDRV  GSEIERMGLV MDRMGSVERM GSGIERMGPL GLDHMASSIE RMGQTMERIG  SGVERMGAGM GFGLERMAAP IDRVGQTIER MGSGVERMGP  AIERMGLSME RMVPAGMGAG LERMGPVMDR MATGLERMGA  NNLERMGLER MGANSLERMG LERMGANSLE RMGPAMGPAL  GAGIERMGLA MGGGGGASFD RAIEMERGNF GGSFAGSFGG  AGGHAPGVAR KACQIFVRNL PFDFTWKMLK DKFNECGHVL YADIKMENGK  SKGCGVVKFE SPEVAERACR MMNGMKLSGR EIDVRIDRNA</p>
<b>Source</b>	E.coli
<b>Target Names</b>	HNRNPM
<b>Protein Names</b>	Recommended name: Heterogeneous nuclear ribonucleoprotein M Short name= hnRNP M
<b>Expression Region</b>	2-730
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	<p>This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid</p>



binding properties. This protein has three repeats of quasi-RRM domains that bind to RNAs. This protein also constitutes a monomer of the N-acetylglucosamine-specific receptor which is postulated to trigger selective recycling of immature GlcNAc-bearing thyroglobulin molecules. Multiple alternatively spliced transcript variants are known for this gene but only two transcripts has been isolated.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.