



recQ Antibody, FITC conjugated

Product Code	CSB-PA17979C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15043
Immunogen	Recombinant Escherichia coli ATP-dependent DNA helicase RecQ protein (2-606AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-dependent DNA helicase RecQ (EC 3.6.4.12), recQ
Immunogen Species	Escherichia coli (strain K12)
Research Area	Microbiology
Target Names	recQ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.