




recQ Antibody

Product Code	CSB-PA17979A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15043
Immunogen	Recombinant Escherichia coli ATP-dependent DNA helicase RecQ protein (2-606AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-dependent DNA helicase RecQ (EC 3.6.4.12), recQ
Immunogen Species	Escherichia coli
Research Area	Microbiology
Target Names	recQ

Image



Western blot
 All lanes: ATP-dependent DNA helicase RecQ antibody at 2µg/ml + DH5a whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 68 kDa
 Observed band size: 68 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.