



SAA1 Antibody, FITC conjugated

Product Code	CSB-PA14047C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DJI8
Immunogen	Recombinant Human Serum amyloid A-1 protein (20-122AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum amyloid A-1 protein (SAA) [Cleaved into: Amyloid protein A (Amyloid fibril protein AA); Serum amyloid protein A(2-104); Serum amyloid protein A(3-104); Serum amyloid protein A(2-103); Serum amyloid protein A(2-102); Serum amyloid protein A(4-101)], SAA1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SAA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.