



SOD1 Antibody, Biotin conjugated

Product Code	CSB-PA02864D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00441
Immunogen	Recombinant Human Superoxide dismutase [Cu-Zn] protein (2-154AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Superoxide dismutase [Cu-Zn] (EC 1.15.1.1) (Superoxide dismutase 1) (hSod1), SOD1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SOD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.