



Recombinant Human Zinc finger protein 133 (ZNF133)

Product Code	CSB-MP026532HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P52736
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAFRDVAVDF TQDEWRLLSQ AQRTLYREVM LENYSNLVSL GISFSKPELI TQLEQGKETW REEKKCSPAT CPADPEPELY LDPFCPPGFS SQKFPMQHVL CNHPPWIFTC LCAEGNIQPG DPGPGDQEKQ QQASEGRPWS DQAEGPEGEG AMPLFGRTKK RTLGAFSRPP QRQPVSSRNG LRGVELEASP AQSGNPEETD KLLKRIEVLG FGTVNCGECG LFSKMTNLL SHQRIHSGEK PYVCGVCEKG FSLKKSARH QKAHSGEKPI VCRECGRGN RKSTLIIHER THSGEKPYMC SECGRGFSQK SNLIHQRTS SGEKPYVCRE CGKGFSQKSA VVRHQRTSLE EKTIVCSDCG LGFSDRSNLI SHQRTHSGEK PYACKECGRC FRQRTTLVNH QRTHSKEKPY VCGVCGHSFS QNSTLISHRR THTGEKPYVC GVCGRGFSK SHLNRHQNIH SGEKPIVCKD CGRGFSQQSN LIRHQRTS EKPMVCGECG RGFSQKSNLV AHQRTHSGER PYVCRECGRG FSHQAGLIRH KRKHSREKPY MCRQCGLGFG NKSALITHKR AHSEEKPCVC RECGQGFLQK SHLTLHQMTH TGEKPYVCKT CGRGFSLKSH LSRHRKTTSV HHRLPVQPDP EPCAGQPSDS LYSL
Source	Mammalian cell
Target Names	ZNF133
Protein Names	Recommended name: Zinc finger protein 133 Alternative name(s): Zinc finger protein 150
Expression Region	1-654
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.