



Recombinant Sorangium cellulosum Lipoyl synthase (lipA)

Product Code	CSB-EP012927SUD-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A9FD62
Product Type	Recombinant Protein
Immunogen Species	Sorangium cellulosum (strain So ce56) (Polyangium cellulosum (strain So ce56))
Purity	≥85% (SDS-PAGE)
Sequence	MAQFSPKPEW LKVRAPGGDT YHHLKETFRK LDLHTVCEE RCPNVGECWR EGTATVMLLG DVCTRGCRFC AVTTGDPRGA VDVREPEHVA RAIARLSLQY VVMTMVNRDD LLDGGAEHVA RTVSRLHALR PDLLIETLVG DFQGHMSAVD MVVDAGPDVF AHNVEVVRRI TRVIRDVRSS YDQSLAVLRR AKERQRRLLA DAAEAGAPAP RRLTKSSIMV GIGETDDEV EALRDLREAG VDIVTIGQYL RPSSKHAPVQ RFVEPETFAA FERAALMGF LYAASAPLVR SSKAAEVFV RSLMDRGGAA LPASPGAAAV EALLEERLAV ARREARLTA ELDPDEPRPP VAPAPASASP ARLVPAASLI RR
Source	E.coli
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.