



# WFDC2 Monoclonal Antibody

<b>Product Code</b>	CSB-MA097901A0m
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	Recombinant Human WAP four-disulfide core domain protein 2
<b>Raised In</b>	mouse
<b>Species Reactivity</b>	Human
<b>Specificity</b>	specific for human WAP four-disulfide core domain protein 2 denatured forms
<b>Tested Applications</b>	ELISA
<b>Form</b>	liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%,Protein G purified
<b>Isotype</b>	IgG1
<b>Clonality</b>	monoclonal
<b>Alias</b>	Epididymal secretory protein E4, Major epididymis-specific protein E4, Putative protease inhibitor WAP5
<b>Product Type</b>	Monoclonal Antibody
<b>Accession NO.</b>	6A1B2