



WFDC2 Monoclonal Antibody

Product Code	CSB-MA097901A0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant Human WAP four-disulfide core domain protein 2
Raised In	mouse
Species Reactivity	Human
Specificity	specific for human WAP four-disulfide core domain protein 2 denatured forms
Tested Applications	ELISA
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%,Protein G purified
Isotype	IgG1
Clonality	monoclonal
Alias	Epididymal secretory protein E4,Major epididymis-specific protein E4,Putative protease inhibitor WAP5
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	6A1B2