



estradiol Monoclonal Antibody

Product Code	CSB-MA004811I0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	estradiol acid-BSA conjugate
Raised In	mouse
Specificity	No significant cross-reactivity or interference was observed
Tested Applications	ELISA
Form	liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%,protein G purified
Isotype	IgG2b
Clonality	monoclonal
Alias	Beta-estradiol, 17beta-estradiol, estrace, oestradiol, dihydrofolliculin, progynon, dihydroxyestrin, dihydrotheelin, diogynets
Product Type	Small Molecules Antibodies
Target Names	E2
Accession NO.	5A9E9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.