





MGMT Antibody

Product Code	CSB-PA11859A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16455
Immunogen	Recombinant Human Methylated-DNA--protein-cysteine methyltransferase protein (4-207AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:10-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methylated-DNA--protein-cysteine methyltransferase (EC 2.1.1.63) (6-O-methylguanine-DNA methyltransferase) (MGMT) (O-6-methylguanine-DNA-alkyltransferase), MGMT
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MGMT

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA11859A0Rb at dilution of 1:20



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA11859A0Rb at dilution of 1:20

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.