



NTH1 Antibody

Product Code	CSB-PA18159A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32356
Immunogen	Recombinant Saccharomyces cerevisiae Neutral trehalase protein (546-751AA)
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neutral trehalase (EC 3.2.1.28) (Alpha,alpha-trehalase) (Alpha,alpha-trehalose glucohydrolase), NTH1, NTH
Immunogen Species	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Research Area	Others
Target Names	NTH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.