



TEP1 Antibody, Biotin conjugated

Product Code	CSB-PA18049D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99973
Immunogen	Recombinant Human Telomerase protein component 1 protein (2-226AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Telomerase protein component 1 (Telomerase-associated protein 1) (Telomerase protein 1) (p240) (p80 telomerase homolog), TEP1, TLP1, TP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TEP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.