



# Recombinant Human Geranylgeranyl transferase type-2 subunit beta (RABGGTB)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-BP019235HU  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | P53611  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | GTPQKDVII KSDAPDTLLL EKHADYIASY GSKKDDYEYC MSEYLRMSGI<br>YWGLTVMDLM GQLHRMNREE ILAFIKSCQH ECGGISASIG HDPHLLYTLS<br>AVQILTLYDS INVIDVNKVV EYVKGLQKED GSFAGDIWGE IDTRFSFCAV<br>ATLALLGKLD AINVEKAIEF VLSCMNFDDG FGCRPGSESH AGQIYCCTGF<br>LAITSQLHQV NSDLLGWWLC ERQLPSGGLN GRPEKLPDVC YSWWVLASLK<br>IIGRLHWIDR EKLRNFILAC QDEETGGFAD RPGDMVDPFH TLFGIAGLSL<br>LGEEQIKPVN PVFCMPEEVL QRVNVQPELV S               |
| <b>Source</b>            | Baculovirus   |
| <b>Target Names</b>      | RABGGTB   |
| <b>Protein Names</b>     | Recommended name: Geranylgeranyl transferase type-2 subunit beta EC=2.5.1.60<br>Alternative name(s): Geranylgeranyl transferase type II subunit beta<br>Short name= GGTase-II-beta Rab geranyl-geranyltransferase subunit beta<br>Short n   |
| <b>Expression Region</b> | 2-331   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |