



Recombinant Danio rerio Methylthioribose-1-phosphate isomerase (mri1)

Product Code	CSB-MP014799DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A9JRE2
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MTLEAIRYRS GSLQILNQLL LPRETVYDEI RSVRDGYEAI KSMKVRGAPA IAIVGCLSLA VELRAGAGAE DLVSFVRDSL CHLTSARPTA VNMGRAAREL MEFTENESME KNTEQLRDSV IGWIEEMLER DVNDNKKIGN YGAQHILSGV PRDSVTILTH CNTGSLATAG YGTALGVVRS LHMLGRLKRL YCTETRPYNQ GARLTAYEAV AEGFPATLIT DSMAALAMRE KSITAVVGA DRVVANGDTA NKVGTYQLAI AAKHHGIPFY VAAPSTSCDL SLESGRDIVI EERPAEELTS INGVPVAAPG IDVWNPAFDV TPHQLITGGI ITELGVFLPS ELQAALTGRL TAL
Source	Mammalian cell
Target Names	mri1
Protein Names	Recommended name: Methylthioribose-1-phosphate isomerase Short name= M1Pi Short name= MTR-1-P isomerase EC= 5.3.1.23 Alternative name(s): S-methyl-5-thioribose-1-phosphate isomerase Translation initiation factor eIF-2B sub
Expression Region	1-353
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.